



# ACCU-View Quick Help



By  
ToupTek Photonics

## ACCU-View Contents

1	ACCU-View .....	1
1.1	Operating System .....	1
1.2	Supported Language .....	1
2	Start to Use ACCU-View.....	2
2.1	Start ACCU-View .....	2
2.2	Open the Camera .....	2
2.3	Adjust Settings of the Camera .....	3
2.3.1	Change the Resolution Setting .....	3
2.3.2	Change the Exposure Time and Gain Setting.....	3
2.3.3	Change the White Balance Setting .....	4
2.3.4	Change the Color Setting .....	4
2.3.5	Change the Frame Rate Setting .....	5
2.3.6	Change the Color mode Setting.....	5
2.3.7	Change the Flip Setting .....	5
2.3.8	Change the Sample Mod Setting .....	6
2.3.9	Change the Power Frequency Setting .....	6
2.3.10	Change the Histogram Setting.....	6
2.3.11	Save & Load Parameter Setting .....	7
2.4	Image Capture and Save .....	7
2.4.1	Single Frame Capture and Save .....	7
2.4.2	Sequence Capture and Batch Save .....	8
2.5	Video Record .....	9
2.6	Open the Browse Window .....	10
2.7	Close Windows .....	11
2.7.1	Close single Windows .....	11
2.7.2	Close All .....	11
3	ACCU-View Windows GUI .....	12
3.1	Function GUI .....	12
3.2	ACCU-View Toolbar .....	13
4	ACCU-View Practical Tool Kit.....	14
4.1	Luminance Calibration Tool .....	14
4.2	Video or Image Calibration Tool .....	14
4.3	Video Overlay Tool .....	15

4.4	Video or Image Measurement Tool .....	15
4.5	Extending Depth of Field (EDoF) Tool .....	16
4.6	Video or Image Stitching Tool .....	16
4.7	ROI Image Select Tool .....	17
4.8	Video or Image Grid Tool .....	17
4.9	Video Watermark Tool .....	18
4.10	Video Marker Tool .....	18
4.11	Line Profile Tool .....	19
4.12	Color Composition Tool .....	19
5	ACCU-View Image Process Modules .....	21
5.1	Filter Shift+F .....	21
5.1.1	Filters: Image Enhance .....	21
5.1.2	Filter: Edge Enhance .....	23
5.1.3	Filters: Morphological .....	24
5.1.4	Filter: Kernel .....	27
5.2	Range••• Shift+R .....	28
5.3	Segmentation••• Shift+S .....	28
5.4	Binary••• Shift+B .....	29
5.5	Emboss••• Shift+E .....	29
5.6	Pseudo Color••• .....	30
5.7	Surface Plot••• .....	30
5.8	Line Profile••• .....	31
5.9	Diffuse•••Shift+D .....	32
5.10	Granuate•••Shift+G .....	32
6	Copyright & Its Term for USage .....	33
6.1	Copyright Statement .....	33
6.2	Terms for usage .....	33

# 1 ACCU-View



ACCU-View is one of the **TOUPTEK PHOTONICS**'s most famous camera control software. It provides functions to fully control the camera and present the videostream processed by Ultra Fine™ color engine at high speed, which includes dedicated pipeline to process the raw data into a realistic scene. Besides, diversified useful tools are provided for specific purpose, such as luminance calibration, measurement, image stitching, extending depth of field, video watermark attachment, color composition, imaging processing and so on. Multi-language mechanism is also realized to support random language, which includes but not limited to English, Chinese, Russian, Turkish, Korea, Polish and so on. Now ACCU-View is widely used in the field of medical microscopic imaging, industrial detection, machine vision, astronomical observation, etc.

ACCU-View is totally compatible with ToupCam full series of digital cameras. With authorized license, ACCU-View software could be used with other cameras, which support Twain or DirectShow interface. ACCU-View is one of the best software in the camera industry, and the United States education department is strongly recommend

## 1.1 Operating System

Microsoft Windows:

- 32bit Windows XP, Vista, 2008, Win7, Win8
- 64bit Windows XP, Vista, 2008, Win7, Win8

Mac:

- OS X 10.6, OS X 10.7, OS X 10.8

Linux:

- Linux 2.6 or above



## 1.2 Supported Language

Standard Language Package:


- 1. Simplified Chinese, 2. Traditional Chinese, 3. English

Optional Language Package:

- 4. German, 5. Japanese, 6. Russian, 7. French, 8. Italian, 9. Polish, 10. Turkish

## 2 Start to Use ACCU-View

### 2.1 Start ACCU-View

1. Install ToupTekACCU-ViewSetup.exe to finish the installation.
2. Double click on the desktop icon  to start ACCU-View;




**Note:** If you are using the Windows 64-bit system, it's recommended to select the x64 ACCU-View to run ACCU-View for maximum performance, even though x86 ACCU-View could also work.

### 2.2 Open the Camera

ACCU-View will detect all of the cameras that your computer has installed (Here, it is **UCMOS03100KPA**, a 3.1M pixel CMOS camera) and will append all the camera names as submenu to the **File>Camera List** menu (Here, the submenu name is "**UCMOS03100KPA**").

Choosing the **File>Camera List>UCMOS03100KPA** will create a video window and start the video stream. The video window will be associated with the name of "**Video [UCMOS03100KPA]**" (i.e., its title bar name will be "**Video [UCMOS03100KPA]**").




There are 4 methods to start the video stream, they are:

File>Camera List	<ul style="list-style-type: none"> <li>◆ Choose the <b>File&gt;Camera List&gt;UCMOS03100KPA</b> (Here, a 3.1M pixel camera is installed) command to create the video window</li> </ul>
Camera Sidebar	<ul style="list-style-type: none"> <li>◆ Click the <b>Camera Sidebar</b> (If it is not activated) and the <b>Camera List</b> to expand the <b>Camera List</b> group (if not expanded). Click the camera name (Here it is <b>UCMOS03100KPA</b>) to create the video window;</li> </ul> 
Toolbar button	<ul style="list-style-type: none"> <li>◆ Click the button 's down arrow on the toolbar to expand the camera list and choose the right camera (Here it is <b>UCMOS03100KPA</b>) to create the video window.</li> </ul> 
Twain Acquire	<ul style="list-style-type: none"> <li>◆ Choose the <b>File&gt;Twain: Select Device...</b> command to select the device <b>UHCCD01400KPA</b> from the <b>Select Source</b> dialog (If never selected before)</li> <li>◆ Choose the <b>File&gt;Twain: Acquire...</b> command. There should be a dialog box like below</li> </ul>


**Note:** ACCU-View supports the connection of one or more cameras to the computer at a time.

## 2.3 Adjust Settings of the Camera

### 2.3.1 Change the Resolution Setting


Control Group	
Function	<ul style="list-style-type: none"> <li>◆ <b>Snap</b> : Click to snap images with specified resolution</li> <li>◆ <b>Record</b> : Click to start to record video</li> <li>◆ <b>Live</b>: Set the <b>Video</b> resolution for preview</li> <li>◆ <b>Snap</b>: Set the <b>Snap</b> resolution for capture</li> </ul>
Remarks	<ul style="list-style-type: none"> <li>◆ <b>Live</b>: Select a resolution in the combobox and videostream will be generated in the specified resolution continuously. Generally a small resolution is selected to obtain faster <b>Video Frame Rate</b></li> <li>◆ <b>Snap</b>: Select a resolution in the combobox and an image in the specified resolution will be obtained when clicking on the Snap button. Generally the maximum resolution is selected for high quality</li> </ul>

### 2.3.2 Change the Exposure Time and Gain Setting


Control Group	
Function	<ul style="list-style-type: none"> <li>◆ When the <b>Exposure &amp; Gain</b> group is expanded, a green rectangle viewfinder marked with <b>Exposure</b> will be overlaid on the <b>Video</b>. This labeled region is taken as the reference region for judging if the average brightness in the region has reached the <b>Exposure Target</b> value. The <b>Exposure ROI</b> rectangle border can be dragged to modify the size</li> <li>◆ Uncheck the <b>Auto Exposure</b> box to switch the the exposure control from <b>Auto Exposure</b> mode to <b>Manual Exposure</b> mode. The <b>Exposure Target</b> slider will be disabled and the <b>Exposure Time</b> slider bar will be enabled for adjustment. Increasing the <b>Exposure Time</b> will increase the brightness and decreasing it will do the opposite. Increasing the <b>Exposure Time</b> has no effect on the frame rate until it reach a specific limit. If the <b>Exposure Time</b> is larger than the limit, the frame rate will be decreased</li> <li>◆ In the <b>Manual Exposure</b> mode at a fixed <b>Exposure Time</b>, the image brightness can also be changed by adjusting the microscope light source</li> <li>◆ Increasing the <b>Gain</b> can also increasing the brightness without having an effect on the frame rate. But <b>Gain</b> with a greater value will also magnify the noise. So only when the microscope light intensity is</li> </ul>

	<p>very weak, drag the <b>Gain</b> slider right to obtain a brighter scene</p>
<p>Remarks</p>	<ul style="list-style-type: none"> <li>◆ The primary user can choose <b>Automatic Exposure</b> mode, through setting the appropriate <b>Exposure Target</b> value, to make the image brightness comfortable to human eye. When the external light source intensity is changed, the software will adjust Exposure Time automatically to make the average brightness in the viewfinder area to the target brightness</li> <li>◆ The other option is that user can uncheck the <b>Auto Exposure</b> and adjust the <b>Exposure Time</b> manually by dragging <b>Exposure Time</b> slider</li> </ul>

### 2.3.3 Change the White Balance Setting


<p>Control Group</p>	
<p>Function</p>	<ul style="list-style-type: none"> <li>◆ Expand the <b>White Balance</b> group and a red rectangular viewfinder marked with <b>White Balance</b> will be overlaid on the video; Drag the viewfinder to a pure white or gray object and click <b>White Balance</b> button to establish the video white balance for future video</li> <li>◆ <b>Temperature</b> and <b>Tint</b> could be used to adjust <b>White Balance</b> setting manually for user defined purpose</li> </ul>
<p>Remarks</p>	<ul style="list-style-type: none"> <li>◆ It's recommended that the temperature value is set to a color temperature value of the light source to obtain a proper white balance effect</li> <li>◆ Users could adjust <b>Temperature</b> and <b>Tint</b> value to obtain a particular effect. For example, set a lower temperature will have a warmer effect on the video</li> </ul>

### 2.3.4 Change the Color Setting

<p>Control Group</p>	
<p>Function</p>	<ul style="list-style-type: none"> <li>◆ <b>Hue</b>: Adjust the Hue value of the video</li> <li>◆ <b>Saturation</b>: Adjust the saturation value of the video</li> <li>◆ <b>Brightness</b>: Adjust the brightness value of the video</li> </ul>

	<ul style="list-style-type: none"> <li>◆ <b>Contrast</b>: Adjust the contrast value of the video</li> <li>◆ <b>Gamma</b>: Adjust the gamma value of the video</li> <li>◆ <b>Defaults</b>: Restore all the color settings to default one</li> </ul>
Remarks	<ul style="list-style-type: none"> <li>◆ <b>Contrast</b> and <b>Gamma</b> adjustment could bring sharper or brighter effect in some condition if users prefer.</li> </ul>


### 2.3.5 Change the Frame Rate Setting

Control Group	
Function	<ul style="list-style-type: none"> <li>◆ This group can be used to adjust the <b>Video Frame Rate</b>. Drag the slider to the right to ensure the fastest <b>Video Frame Rate</b> if your computer USB hub supports</li> </ul>
Remarks	<ul style="list-style-type: none"> <li>◆ If <b>Video</b> window appears black and fails to display, drag the slider to the left one step a time until the video window could display video normally. This condition happened when the USB transportation capability of hardware on PC side is not strong enough</li> </ul>

### 2.3.6 Change the Color mode Setting


Control Group	
Function	<ul style="list-style-type: none"> <li>◆ <b>Color</b>: Set the video to <b>Color</b> mode</li> <li>◆ <b>Gray</b>: Set the video to <b>Gray</b> mode</li> </ul>
Remarks	<ul style="list-style-type: none"> <li>◆ <b>Gray</b> images can be used in applications without color requirement</li> </ul>

### 2.3.7 Change the Flip Setting


Control Group	
Function	<ul style="list-style-type: none"> <li>◆ <b>Horizontal</b>: If the <b>Video</b> on the screen appears in different horizontal direction, check this to set the video direction to the right one</li> <li>◆ <b>Vertical</b>: If the <b>Video</b> on the screen appears in different vertical direction, check this to set the video direction to the right one</li> </ul>
Remarks	<ul style="list-style-type: none"> <li>◆ If the video has a orientation deviation of 90 or 270 degrees, please rotate the camera relative to the sample to adjust the video orientation</li> </ul>




### 2.3.8 Change the Sample Mod Setting

Control Group	
Function	<ul style="list-style-type: none"> <li>◆ <b>Bin</b>: Pixel binning refers to the method of combining (averaging) adjacent pixels of the same color. Since <b>Bin</b> will take some time, the frame rate in <b>Bin</b> mode is slower than that in <b>Skip</b> mode normally</li> <li>◆ <b>Skip</b>: Also called <b>Decimation</b>, means that a certain amount of pixels is not read out but skipped (horizontally, vertically or in both axes). This reduces resolution of the resulting image to obtain a high frame rate with a side effect of introducing sub sampling artifacts</li> </ul>
Remarks	<ul style="list-style-type: none"> <li>◆ <b>Bin</b>: Good image quality, spend a little more time on calculation while the frame rate will be reduced</li> <li>◆ <b>Skip</b>: Fast frame rate while the image quality is slightly lower</li> </ul>


### 2.3.9 Change the Power Frequency Setting

Control Group	
Function	<ul style="list-style-type: none"> <li>◆ This function can eliminate stroboscopic phenomenon between different countries and areas for different power frequency, please select it according to your national local power characteristics</li> </ul>
Remarks	<ul style="list-style-type: none"> <li>◆ It is <b>50Hz</b> for China</li> <li>◆ If the stroboscopic phenomenon still appears after choosed the right power frequency, please set the <b>Exposure Time</b> to the integer times of <math>10\text{ms}(50\text{Hz})/8.333\text{ms}(60\text{Hz})</math>. This can eliminate the stroboscopic phenomenon completely</li> </ul>

### 2.3.10 Change the Histogram Setting

Control Group	
Function	<ul style="list-style-type: none"> <li>◆ A <b>Histogram</b> illustrates how pixels in an image are distributed in each color channel. A <b>Histogram</b> can help you determine whether an image has enough detail to make a good correction</li> <li>◆ There are two lines on right and left respectively. The pixels will always be destributed between the two lines. Resposition the lines will change the distribution of the pixels</li> </ul>
Remarks	<ul style="list-style-type: none"> <li>◆ Drag the left vertical line to the right will increase the image contrast while drag the right vertical line to left will will increase the image background brightness</li> </ul>

### 2.3.11 Save & Load Parameter Setting


Control Group	
Function	<ul style="list-style-type: none"> <li>◆ After the light source and the camera parameters are adjusted to make the video satisfied, the camera control settings can be saved by clicking the <b>Save</b> button. The saved parameter can be <b>Load</b> for the later use to resume the same scenario</li> </ul>
Remarks	<ul style="list-style-type: none"> <li>◆ The parameters could also be save, load, overwritten</li> </ul>

## 2.4 Image Capture and Save

### 2.4.1 Single Frame Capture and Save

Click the **Snap** button on **Camera Control Sidebar**, an image can be obtained in a specified resolution.

**Note:** Using the ROI tool to select the region that you are interested in, the snapped image will be only the ROI area.

To save the current captured image, click **Save**  button on the toolbar will invoke a dialog called **Save as**. Fill the content in the dialog and click **Save** will save the captured image. ACCU-View supported:

Window Bitmap(\*.bmp, \*.dib, \*.rle)、JPEG(\*.jpg, \*.jpeg, \*.jpe, \*.jif, \*.jfif) Portable Network Graphics(\*.png)、Tag Image File Format(\*.tif, \*.tiff)、Compuserve GIF (\*.gif)、PCX(\*.pcx) Targa(\*.tga) JBIG(\*.jbg)、ACCU-View File Type(\*.tft)

**Note:** For WindowBitmap(\*.bmp, \*.dib, \*.rle), Compuserve GIF (\*.gif), PCX(\*.pcx), Targa(\*.tga), JBIG(\*.jbg), ACCU-View File Type(\*.sft), the **Options** button is not enabled to indicate that these formats do not have optional settings.

For JPEG(\*.jpg, \*.jpeg, \*.jpe, \*.jif, \*.jfif)、Portable Network Graphics(\*.png)、Tag Image File Format(\*.tif, \*.tiff), the **Option** button is enabled. The corresponding dialogs are shown as below(Fig.2~Fig.4) :




Fig.2




Fig.3




Fig.4




Fig.5

Image quality	◆ In case to save an image in JPEG format (*.jpg), Image quality in the edit box is available for adjustment. The values range from 0 to 100. Default value:75
Progressive	◆ The default is unchecked
Optimize coding	◆ The default is unchecked
Smoothing	◆ The values range from 0 to 100. Default value:0
Save these setting as defaults	◆ When saving a file, the current settings will be saved as defaults for the next file save operation if this box is checked
Interlaced	◆ The default is unchecked
Appended pages	◆ Determine whether the current image will be saved in multiple pages style or not
Compressions	◆ TIFF format can be chosed in different compression ways, default value: LZW (default)
Image quality	◆ In case of choosing Compressions as JPEG, the Image quality is available or adjustment. The values range from 0 to 100. Default value: 75
Reset	◆ Set to a specified value
Save these settings as defaults	◆ The current settings will be saved as defaults for the future

### 2.4.2 Sequence Capture and Batch Save

Choose Capture>Time>Lapse menu to set the Time-lapse(Auto Capture) parameter. (See Fig.1)




Fig.1

Sequence capture is also called **Time-Lapse Capture**, This command can capture a sequence of pictures in a predefined interval. User can specify the directory and naming format for the captured images, precisely set the interval (Time slot, from 2 to 3600 seconds) and the total number of images.

When a lot of images are capture, the **File>Batch Save...** will be enabled.

The captured files can be saved at a time by choosing **File>Batch Save...** menu(Fig.5).

## 2.5 Video Record




Fig.6






Fig.7

Clicking the **Record** button  on the **Camera** control sidebar in the **Capture & Resolution** Group will bring up the **Video Format** dialog (Fig.6) . Check the **wmv/asf** or **avi** and click **Next** to invoke the **Video File** dialog as shown in Fig.7.

Fill the **Video File Name** and **Directory**, clicking **Next** will invoke the **Encoder** dialog (Fig.8).

Here you can select the **Encoder** format, set the **Bitrate**(Kbps), **Quality**(1-100) and **Key Frames Spacing** (1-30), clicking **Next** will invoke the **Display Information** dialog(Fig.9)

Here you can enter **Title**, **Author**, **Copyright** and **Description** to their fields. Clicking **Back** to return to the Encoder dialog, or **Next** to the next step (Fig.10).






Fig.9





Fig.10

A dialog called **Start to Record** will be brought up(Fig.10). Here you can check **Time Limit** (from 1 to 1440 minutes), enter recording time (If checked) and **Time-lapse** (1-100). There is a **Summary** that displays all the options. Click **Finish** to start to record the video.

After the **Video Record** is started. The **Record**  button on the **Camera Sidebar** will become . Clicking  will stop the recording process.

## 2.6 Open the Browse Window

Choose the **View>Browse** menu from the **View** menu or click the **Browse** toolbar button  to browse images under the specified directory on the **Folders Sidebar**;

Click the **Folders Sidebar** to activate it and double-clicking the listed directory in the **Folders Sidebar** will create a **Browse** window as shown below(Fig.11).





Fig.11


After creating the **Browse** window, ACCU-View will display the **Browse** window that looks like windows explorer. The child window on the left part of the **Browse** window called **Folders Sidebar** and can be used to

relocate the directory on the hard disk. Images in the current directory can be displayed in **Large Icons** or **Small Icons** mode on the right side of the **Browse** window upon your setting.

When the **Browse** window is active, the **Edit** and **View** menus will be changed to the following styles. There are shown as below.



The **Edit** menu(Fig.12)



The **View** menu(Fig.13)

These two menus can be used for the basic setting of the **Browse** window. For example, image file order can be set in **Forward** or **Reverse** format according to **Sort by Name, Type, Size, Width** or **Height** et al.

## 2.7 Close Windows

### 2.7.1 Close single Windows

Click on the window top-left corner to close the corresponding or double-click the name of the window.


### 2.7.2 Close All

Choose **Windows>Closes All** menu to close all of the pictures opened or captured inside ACCU-View frame.



# 3 ACCU-View Windows GUI

## 3.1 Function GUI



A: ACCU-View; B: Menu; C: ACCU-View toolbar D: Camera Sidebar ;

E: Folders Sidebar ; F: Undo/Redo Sidebar ; G: Layer Sidebar ;

H: Measurement Sidebar ; I: Measurement Sheet; J: Statusbar;

K: Auto Hide button L: Horizontal ruler; M: Vertical ruler

N: Frame Rate O: Frames captured P: Current Video sizes

Q:Selected microscope **Magnification** R:Current **Unit**;

AA: **Sidebar** right mouse button context menu;

AB: **Video** window right mouse button context menu;

AC: **Browse** window right mouse button context menu;

AD: **Image** window right mouse button context menu;

AE: **Frame** window right mouse button context menu;

AF: Double-click bring up **Video Properties** dialog;

AG: Double-click bring up **Magnification** dialog;

AH: Double-click bring up **Measurement** dialog;

AI: **Horizontal Ruler** or **Vertical Ruler** right mouse button context menu

### 3.2 ACCU-View Toolbar






Item	Function	Item	Function
1	Open ( Ctrl+O)	16	Parallel
2	Save( Ctrl+S)	17	Two Parallel
3	Browse(Ctrl+B)	18	Vertical
4	Camera List	19	Rectangle
5	Unit	20	RoundRect
6	Magnification	21	Ellipse
7	Zoom	22	Circle
8	Video Source Properties	23	Annulus
9	Video Stream Format	24	Two Circles
10	Video/Image Select	25	Arc
11	Track	26	Text
12	Object Select	27	Polygon
13	Angle	28	Calibration
14	Point	29	Gray Calibration
15	Line	30	Manual Fusion




# 4 ACCU-View Practical Tool Kit

## 4.1 Luminance Calibration Tool


Control Dialog	
Function	<ul style="list-style-type: none"> <li>◆ Switch the camera from the <b>Auto Exposure</b> mode to the <b>Manual Exposure</b> mode in the <b>Exposure &amp; Resolution</b> group</li> <li>◆ Use ROI <b>Image Select</b> button  to select a ROI on the background</li> <li>◆ Click <b>Gray Calibration</b> button  on the toolbar and the above dialog appeared</li> <li>◆ Adjust light source brightness to make the <b>Average Gray</b> value to the desired value</li> </ul>
Initial Set	<ul style="list-style-type: none"> <li>◆ Recommended <b>Average Gray</b> value: 190~205</li> </ul>

## 4.2 Video or Image Calibration Tool




Function	<ul style="list-style-type: none"> <li>◆ Set the <b>Live Resolution</b> to the maximum resolution in <b>Capture &amp; Resolution</b> group</li> <li>◆ Set the <b>Zoom</b> ratio to 100%</li> <li>◆ Click <b>Calibration</b> button on the toolbar</li> <li>◆ Drag the reticule to align the calibration tails</li> <li>◆ Input the objective lens <b>Magnification</b>, <b>Actual Length</b> of the calibration tails, and then click <b>OK</b> button to end the calibration</li> <li>◆ *The <b>Export</b> button is used to export the calibrations to a file. This will save the calibration parameters and can be <b>Imported</b> again for the future application</li> </ul>
Initial Set	<ul style="list-style-type: none"> <li>◆ *Make sure to use the right <b>Magnification</b> calibration file to match the actual lens' magnification</li> </ul>



### 4.3 Video Overlay Tool

<p>Control Dialog</p>	
<p>Function</p>	<ul style="list-style-type: none"> <li>◆ Choose <b>Setup&gt;Video Overlay</b> menu to overlay the current <b>Magnification</b>, <b>Scale Bar</b>, <b>Date Time</b> and <b>Clarity Factor</b> on the video</li> </ul>
<p>Initial Set</p>	<ul style="list-style-type: none"> <li>◆ *The <b>Overlays</b> can be overlaid on the captured image by checking <b>Capture with Marker</b> and <b>Watermarker</b> box on the <b>Options&gt;Preferences&gt;Capture</b> dialog</li> </ul>



### 4.4 Video or Image Measurement Tool

<p>Control Dialog</p>	
<p>Function</p>	<ul style="list-style-type: none"> <li>◆ <b>Measurement</b> tool includes common two-dimensional plane measurement command. Choose suitable tool according to specific application</li> <li>◆ Click the measure anchor and click the left button of the mouse to finish the measurement. Click the right button of the mouse will abort the current measurement. The <b>Measurement Object</b> can be moved, deleted or modified by selecting it with the <b>Object Select</b> tool</li> <li>◆ Choose the <b>Export to Image</b> in <b>Layer</b> menu to combine the <b>Measurement Objects</b> with the background image to form a new background image</li> <li>◆ All the measure anchors are on a new layer over the background image, this will not damage the original image. After the object is selected, click <b>Delete</b> on keyboard to delete them</li> </ul>
<p>Initial Set</p>	<ul style="list-style-type: none"> <li>◆ The <b>Object</b> properties can be predefined by <b>Options&gt;Measuremen&gt;General</b> or <b>Object</b> page</li> </ul>



## 4.5 Extending Depth of Field (EDoF) Tool

<p>Control Group</p>	
<p>Function</p>	<ul style="list-style-type: none"> <li>◆ Click <b>Manual Fusion</b>  on the toolbar and click <b>Capture</b> to capture the image for Fuse operation</li> <li>◆ Finely and evenly adjust the Z-Axis and capture image once a time for each adjustment. Click <b>Manual Fusion</b> after finish capturing images. Wait for a while, and then a new <b>EDoF</b> image window will be created</li> <li>◆ This <b>EDoF</b> extracts the clearest part from the different focal planes to take the place of the fuzzy part to obtain a clear image. Please pay attention to the changes of every adjustment and make sure that each image has a different clear part. To avoid the misalignment, please don't change the XY direction in the <b>EDoF</b> process.</li> </ul>
<p>Initial Set</p>	<ul style="list-style-type: none"> <li>◆ Make sure that each image should have a different clear part. It is not necessary to capture too many images.</li> <li>◆ Rotate the fine knob evenly along the same direction in the <b>EDoF</b> process</li> <li>◆ Images without clear parts will make the final picture fuzzy</li> </ul>


## 4.6 Video or Image Stitching Tool

<p>Control Dialog</p>	
<p>Function</p>	 <ul style="list-style-type: none"> <li>◆ Uncheck the <b>Auto Exposure</b> box to set the <b>Exposure</b> mode to <b>Manual Exposure</b> mode and make the image clear.</li> </ul>

Initial Set


- ◆ Choose **Process>Stitching** (  ) menu , click **Snap** to capture image and then move the sample forward along a direction, click **Snap** and click capture again to repeat the operations until the end. After all the images have been captured, Click **OK** to begin image stitching process and after a while, a stitched image will be generated in a new window
- ◆ Make sure that the image brightness will not change significantly during the movement
- ◆ \*Each moving distance should not exceed 75% of window content, which means that there should have 25% overlap region between every 2 adjacent images.
- ◆ Click **Image Select** button  to select the desired area and click **Image>Crop** (Shortcut Key: Shift+C) to remove the ragged black edge (caused by the movement that is not in an extremely straight line). Finally, choose save to finish the image **Stitching** process.

## 4.7 ROI Image Select Tool



Function



Initial Set

- ◆ **Image Select**  button on the toolbar has many applications.
- ◆ Select a region in the video window and click **Snap** to obtain a selected **ROI** image
- ◆ Select a region in the video window and click **Gray Calibration** tool to calibrate the video gray to a desired value for the **Luminance** calibration.
- ◆ Select a region in the video window will enable the **Copy** button
- ◆ \*Click **Image Select** button to activate it and click again to close it

## 4.8 Video or Image Grid Tool


Control

Dialog




Function	<ul style="list-style-type: none"> <li>◆ Operating method: Choose <b>View&gt;Ruler and Grid&gt;Grids&gt;Auto Grids</b> menu</li> <li>◆ It can be used as the virtual cell counting chamber.</li> <li>◆ *<b>Auto Grids</b>: The grids density is changed according to the image resolution</li> </ul>
Initial Set	<ul style="list-style-type: none"> <li>◆ *<b>Manual Grids</b>: Drag the horizontal or vertical arrow at the top left of the window to a desired region manually will overlay the a grid line on the image.</li> </ul>


## 4.9 Video Watermark Tool

Control Dialog	
Function	<ul style="list-style-type: none"> <li>◆ See the detailed information in <b>Setup&gt;Watermark</b> section in the ACCU-View help file</li> </ul>
Initial Set	<ul style="list-style-type: none"> <li>◆ The watermark file must be in 24bits BMP</li> <li>◆ To intensify the <b>Watermark</b> effect, the <b>Watermark</b> image could perform binary and invert processing. For example, if a ruler image is taken as the <b>Watermark</b>, you should binarize the ruler first and then invert the calibration tails, as a result, the black tails changed into white and the white background changed into black.</li> </ul>


## 4.10 Video Marker Tool




Control Dialog	
Function	<ul style="list-style-type: none"> <li>◆ Operating method: Choose <b>Setup&gt;Video Overlay</b> menu and click the <b>Marker</b> page</li> </ul>
Initial Set	<ul style="list-style-type: none"> <li>◆ * The <b>Marker</b> overlaid on the video window can be took as real <b>Cross</b> in the eyepiece</li> </ul>

## 4.11 Line Profile Tool

<p>Control Dialog</p>	
<p>Function</p>	<ul style="list-style-type: none"> <li>◆ Operating method: Choose <b>Process&gt;Line Profile</b> menu</li> </ul>
<p>Initial Set</p>	<ul style="list-style-type: none"> <li>◆ *Click the <b>Line</b> button on the toolbar to draw a <b>Line</b>, choose <b>Process&gt;Line Profile</b> menu will invoke the <b>Line Profile</b> dialog. The R,G and B value under the line will be plotted in the <b>Image Line Profile</b> window. It reflects the linear brightness color value distribution under the <b>Line</b> at real time. The <b>Profile</b> line will be changed according to the <b>Line</b> position or when the <b>Line</b> state</li> </ul>

## 4.12 Color Composition Tool

<p>Control Dialog</p>	
<p>Function</p>	<ul style="list-style-type: none"> <li>◆ Operating method: Choose <b>Process&gt;Color Composite</b> menu</li> <li>◆ When and only when a gray image is available, the <b>Color Composite</b> menu will be enabled, so pls change the captured image to gray image first.</li> </ul>
<p>Initial Set</p>	<ul style="list-style-type: none"> <li>◆ *Choose the right color in the dye list according to the dye name and bring the grey image back to colorful. If the dye list do not have the required dye, please define it with the right color first</li> </ul>

<p>The Original Image</p>		
<p>Colorized by Pseudo Color</p>		
<p>Color Composed Image</p>	<p>◆ During the <b>Color Composite</b> process, select the images, adjust its brightness (B), contrast ratio(C) and gamma value (G) separately. Make sure that the size and content of the images fit well.</p>	



# 5 ACCU-View Image Process Modules

## 5.1 Filter Shift+F


Operating Method: Choose **Process>Filter** Menu.

TouView provides various kinds of filters through which you can acquire ideal images efficiently. Generally, the filters are divided into two categories: **Convolution** filters and **Non-convolution (Morphological)** filters.

No matter what categories the filter is, different kernels take different effect on the filtered results despite the same filter category. ACCU-View also allows access to custom filter kernel so that satisfactory results can be achieved during image processing. The filter dialog involves four property sheets or tabs, each of which has its own filter. Either one you choose, you can find real-time effect in preview window. Click **OK** or **Apply**, filtered results are always written to the active image. The difference for these two buttons is, while writing results, click **OK**, dialog closed, click **Apply**, dialog remains. For the results that already written to images, **Edit>Undo** command will cancel the applied operations if you are not satisfied.

### 5.1.1 Filters: Image Enhance

Shown as below, this property page includes filters as follows:



Item	Description
Low Pass	◆ The <b>Low Pass</b> filter replaces the center pixel with the mean value in its neighborhood. Check this filter can soften image and remove noise by eliminating high-frequency information (this has the effect of blurring sharp edges)
High Pass	◆ The <b>High Pass</b> filter enhances high-frequency information by replacing the center pixel with a convolved value calculated from a certain kernel
Gauss	◆ Similar to <b>Low Pass</b> filters, <b>Gauss</b> filter can soften image by eliminating high-frequency information, but not as effective as <b>Low Pass</b> . It softens image with its Gaussian formula



High Gauss	◆ Similar to <b>Sharpen</b> filter, <b>High Gauss</b> enhances detailed information of images, but with less noise and a Gaussian curve type of kernel
Equalization	◆ The <b>Equalization</b> filter enhances image contrast based on the histogram algorithm
Sharpness	◆ The <b>Sharpness</b> filter enhances detailed information of images
Median	◆ <b>Median</b> filter removes impulse noise from an image. It first sorts out the neighborhood pixels according to luminance value, then replaces the current pixel with the <b>Median</b> value
Rank	◆ Similar to <b>Median</b> filter, <b>Rank</b> filter also removes impulse noise from an image. The pixels in the kernel are ranked by order of intensity, and the pixel in that range at the rank percentage is chosen for comparison. For example, in a 5x5 kernel, there are 25 pixels. A rank percentage of 95% would choose second-brightest pixel for comparison. If the difference between the selected pixel and the center pixel is greater than the threshold value, the Rank filter replaces the value of the center pixel with the value of the selected pixel

The functions of filters are in relation with the factors below:

To **Image Enhance** filters, the **shape** and **size** of kernel is of great importance;

To **Equalization** filters, their options relate to histogram **Equalizations**, among them is **Local Histogram Equalization** which modifies the contrast of an image based on the pixel values in a small window surrounding each pixel


Item	Description
3x3	◆ Check <b>3 x 3</b> kernel will produce a more subtle filtering effect
5x5	◆ Check <b>5 x 5</b> kernel will produce a moderate filtering effect
7x7	◆ Check <b>7 x 7</b> kernel will produces a more extreme filtering effect
Passes	◆ Filter applied times on the image. Each operation is based on the results of previous operation. So more times applied, more amplified the effect will be
Strength	◆ Applied value reflecting how much of the filtering effect on the image. The range is from 1 to 10. Value 10 specifies the full strength (100%) of the filtered result applied to each pixel. Value 1 indicates that only 10% of the difference between the filtered pixel value and the original pixel value should be applied
Rank	<p>◆ This value specifies which pixel in the sorted array will be used to replace the center pixel. Pixels in the array will be sorted in ascending order. The pixels are indexed from 0 to <b>Kernel Size x Kernel Size-1</b>. In the pixel index 0 corresponds to the lowest pixel value</p> <p>◆ The <b>Rank</b> will be specified in terms of a percentage of the indexes (<b>Kernel Size x Kernel Size-1</b>). 0% rank means the lowest index (lowest gray value), 50% <b>Rank</b> means the middle of the array. and 100% rank means the highest index (highest gray value)</p>

The **Equalization** filter enhances image contrast based on the histogram algorithm, it has the following options.

Item	Description
Local Histogram Equalization	◆ <b>Local Histogram Equalization</b> equalizes an entire image or partial image according to the size of window
Window	◆ Image pixels statistics (min, max, histogram, mean, standard deviation, etc.) will be calculated on a small <b>Window</b> of the image. These measurements are then used to derive the local contrast for that area of the image. In short, an area of <b>Window x Window</b> around each pixel is all that is considered when modifying the intensities in the image. Larger Window produces smoother results, while small Window track small details more closely
Best Fit	◆ Choose <b>Best Fit</b> command to optimize the values for the particular image. The results are achieved by stretching the local histogram to maximize the contrast between the brightest and darkest pixels in the local window region
Linear	◆ This option distributes the histogram linearly across the intensity scale. This function produces a high contrast image with the highest possible dynamic range
Logarithmic	◆ This option concentrates the histogram at the low end of the scale. This function produces a high contrast image with little dynamic image. It will tend to darken the image overall. It is useful for increasing the contrast in a very light image
Exponential	◆ This option concentrates the histogram at the high end of the scale. This function produces a high contrast image with little dynamic image. It will tend to lighten the image overall. It is useful for increasing the contrast in a very dark image

### 5.1.2 Filter: Edge Enhance

Shown as below, this property page includes filters as follows:



Item	Description
Sobel	◆ Check this filter to enhance the principal edges in an image. (The Sobel applies a mathematical formula to a 3x3 neighborhood to locate and highlight its edges)
Roberts	◆ Check this filter to enhance fine edges in an image. The Roberts filter is not a convolution filter. It applies a mathematical formula upon a 4 x 4 neighborhood to produce its effect. The upper left pixel in the neighborhood is the one that is replaced
Sculpt	◆ Check this filter to apply a sculpted effect on the image
Laplace	◆ Select this filter if you want to enhance all the edges in an image
Variance	◆ Select this filter if you want to detect and emphasize edges and textures. The Variance filter substitutes a pixel with the standard deviation for its neighborhood
Horizontal	◆ Check this filter to detect and emphasize horizontal edges
Vertical	◆ Check this filter to detect and emphasize vertical edges

If one of the Edge filters (Laplace, Variance, Horizontal and Vertical Edge filters) has been checked, the options will relate to kernel size, passes, and filtering strength. The following options will be displayed


Item	Description
3 x 3	◆ Check 3x3 kernels to produce a more subtle filtering effect
5 x 5	◆ Check 5x5 kernels to produce a moderate filtering effect
7 x 7	◆ Check 7x7 kernels to produce a more extreme filtering effect
Passes	◆ Enter the number of times that the filter will be applied to the image. When a filter is applied multiple times, its effect is amplified by each pass. An image that has been softened by one pass of the Image Enhancement Filter, will be softened further by a second pass
Strength	◆ Enter a value from 1-10 that reflects how much of the filtering effect to apply to the image. A value of 10 specifies that the full strength (100%) of the filtered result will be applied to each pixel. Values less than 10 cut the full weight of the filter - a value of 1 indicates that only 10% of the difference between the filtered pixel value and the original pixel value should be applied, a value of 2 indicates that 20% of the difference should be applied, and so forth

**Note:** If you have selected Sobel, Phase or Roberts, no options are available.

If you have selected Variance, the Passes and Strength options are not available

### 5.1.3 Filters: Morphological

Shown as below, this property page includes filters as follows:



Item	Description
Erode	◆ Check this morphological filter if one wants to modify the size of objects in the image. The <b>Erode</b> filter erodes the edges of bright objects and enlarges the edges of dark ones
Dilate	◆ Check this morphological filter if one wants to modify the size of objects in the image. The <b>Dilation</b> filter dilates bright objects and erodes dark ones
Open	◆ Check this morphological filter if one wants to modify the shape of objects in the image. Assuming the image contains bright objects on a dark field, the <b>Open</b> filter will smooth object contours, separate narrowly connected objects, and remove small dark holes
Close	◆ Check this morphological filter if one wants to modify the shape of the objects in the image. Assuming the image contains bright objects on a dark field; the <b>Close</b> filter will fill gaps and enlarge protrusions to connect objects that are close together
Tophat	◆ Check this filter to detect and emphasize points, or grains, that are brighter than the background. There are 3 kernel sizes for this processing. Click the radio button to change the kernel size to the value that most closely matches the size of the grains to detect
Well	◆ Check this filter to detect and emphasize points, or grains, that are darker than the background. There are 3 kernel sizes for this processing. Click the radio button to change the kernel size to the value that most closely matches the size of the grains to detect
Gradient	◆ Check this filter to enhance edges in an image
Watershed	◆ Check this filter to separate objects that are touching. The <b>Watershed</b> filter erodes objects until they disappear, then dilates them again, but will not allow them to touch. The <b>Watershed</b> filter will not operate upon <b>True Color</b> images. If one wants to separate objects in a <b>True Color</b> image, he must first convert it to <b>Gray Scale</b> (see : <b>Image&gt;Gray Scale</b> )
Thinning	◆ Check this filter to reduce an image to its skeleton. When choosing this filter, one must set the threshold that determines whether a pixel is part of the subject, or part of the background (see <b>Options</b>

	below). The Thinning filter will not operate upon <b>True Color</b> images. If one wants to thin a <b>True Color</b> image, he must first convert it to <b>Gray Scale</b>
Distance	<ul style="list-style-type: none"> <li>◆ The <b>Distance</b> filter is used to show the distances of pixels within blobs to the outer boundaries of those blobs. After applying the distance filter, the background will be black (i.e. pixels with value 0). Only the area within the blobs will have non-zero values (will be white). The values of each pixel within the blob will be a count of the shortest distance from that pixel to the edge of the blob. Thus, all pixels along the blob's border will have a value of 1 (since they are one pixel away from the edge of the blob); pixels that are a distance of 2 from the border will have the value 2, and so on. This creates a distance map of the image. The Distance filter will not operate upon <b>True Color</b> images. If one wants to use the <b>Distance</b> filter with a <b>True Color</b> image, he must first convert it to <b>Gray Scale</b></li> </ul>

If **Erode**, **Dilate**, **Open**, or **Close** filters is checked, the options will relate to the kernel size and shape. The following options will be presented:

Item	Description
2 x 2	◆ Check to use the 2x2 square kernel configurations
3 x 1 Row	◆ Check to use the 3x1 row kernel configuration
1 x 3 Column	◆ Check to use the 1x3 column kernel configuration
3 x 3 Cross	◆ Check to use the 3x3 cross kernel configuration
5 x 5 Circle	◆ Check to use the 5x5 circular kernel configurations
7 x 7 Circle	◆ Check to use the 7x7 circular kernel configurations. This is a two-pass filter, accomplished using a 5 x 5 circle followed by a 3x3 cross
11 x 11 Circle	◆ Check to use the 11 x 11 circular kernel configurations. This is a three-pass filter, accomplished using a 5 x 5 circle followed by another 5 x 5 circle, followed by a 3 x 3 cross
Passes	◆ Set the number of times iterate the filter

If the **Tophat**, **Well**, or **Gradient** filter is selected, the options will relate to kernel size and shape. The following options will be presented:

Item	Description
3x3	◆ Check to use the 3x3 square kernel configurations
5x5	◆ Check to use the 5x5 square kernel configurations
7x7	◆ Check to use the 7x7 square kernel configurations


If **Distance**, **Thinning**, or **Watershed** filter is checked, the options will relate to the threshold. The following option will be presented:

Item	Description
Threshold	Enter a percentage value from 1-100 that specifies the intensity value to <b>Threshold</b> binarize the

image. For example, a **Threshold** of 50% on a **Gray Scale** image would set all values  $\leq 127$  to 0 (black) and all values  $\geq 128$  to the maximum value for that image class (white)

### 5.1.4 Filter: Kernel

The **Kernel** page allows editing kernel files for both **Convolution** and **Morphological** filters. The dialog shows as below:



**Note:** The **HiPass**, **LoPass**, **Laplace** and **Unsharp** kernel files are used by the **HiPass**, **LoPass**, **Laplacian** and **Sharpen** options listed in the **Filter** window (i.e., there is no difference between selecting one of these kernel files and selecting its option button in the **Filter** window -- the two methods ultimately do the same thing). Because these kernel files are essential to the operation of these filtering options, they must not be deleted or renamed.


Item	Description
Filter type	<ul style="list-style-type: none"> <li>Check to modify the kernel for a selected <b>Filter type</b>, either <b>Convolution</b> or <b>Morphological</b> filters</li> </ul>
Edit	<ul style="list-style-type: none"> <li><b>Name</b> :This list box contains the name of the selected kernel file. If one wants to save the modified kernel file to the same file,leave it as it is. If one wants to save the file to a new location, enter the new filename here</li> </ul>
	<ul style="list-style-type: none"> <li><b>Kernel Size</b>: Click the spin buttons or enter the number to change the size of the kernel. Either direction may take into account one to nine pixels. As one modifies the <b>Kernel Size</b>, the shape of the kernel representation changes accordingly. In the center of the dialog, there are white boxes containing coefficients that will be multiplied with each pixel that will be taken into account by the filter kernel. One can change any coefficient by clicking on it and adjust it as desired</li> </ul>
	<ul style="list-style-type: none"> <li><b>Fill</b>: Click this button to fill every element of the kernel with a particular value. The <b>Fill</b> kernel dialog appears. One may enter a value between 0 and 10. Using the <b>Fill</b> button is useful for setting all coefficients to the same value. One may then change the coefficients that require a different value</li> </ul>
	<ul style="list-style-type: none"> <li><b>Offset</b>: The pixel whose value is being modified is usually the center-most pixel. One may, however, designate any pixel. ACCU-View signals the pixel to be changed by putting a box around it. Choose X</li> </ul>

	and Y <b>Offset</b> spin buttons to apply
New...	◆ Click to create a new filter kernel. The <b>Edit Kernel</b> dialog will appear. The functions of the dialog are the same way as the dialog for <b>Edit...</b> described above), with the exception that the file name for the new kernel file must be provided
Delete	◆ Click to <b>Delete</b> the selected filter kernel file

## 5.2 Range... Shift+R

Operating method: Choose **Process>Range** menu

The **Range** command allows set the intensity levels of the image to increase the contrast and enhance the display in low-light situations.




Item	Description
High & Low end	◆ The <b>Range</b> dialog shows the Histogram of current image on which two vertical marks stand for high and low end of strength. The cursor can move the marks. While looking at a colorful image, the Histogram will show red, green and blue colors through red line, green line, blue line respectively
2 Edit Controls	◆ The edit controls specify values of strength extreme. Click the up or down arrow on the spin button can increase or decrease the value. Those values that from 0 to low end are black and the value that ranges from high end to maximum of scale is white
Reset	◆ The <b>Reset</b> button allows to reset the black and white levels to the high and low ends of the dynamic <b>Range</b> . <b>Reset</b> only has effect on the displayed range, <b>Reset</b> all functions will display original settings
Best Fit	◆ The <b>Best Fit</b> button automatically sets the intensity levels to the <b>Best Fit</b> . It instructs ACCU-View to optimize the brightness and contrast values for the particular image
Invert	◆ The <b>Invert</b> button reverses the color of the image
Update	◆ <b>Update</b> will refresh the display <b>Histogram</b> of the current dialog

## 5.3 Segmentation... Shift+S

Operating method: Choose **Process>Segmentation** menu.

**Segmentation** process is a Histogram based models for identifying certain colors in the segmented range. The segmented range can either be selected by two vertical bars with mouse, or directly editing value through the two


edit controls. The colors used are specified by the system palette. For the 24 bit true color image, the Segmentation operation can be performed in either Red(R), Green(G), or Blue(B) channel separately.



### 5.4 Binary... Shift+B

Operating method: Choose **Process>Binary** menu.

**Binary** is a kind of gray level process. If the gray of the pixel is greater than the given threshold, the pixel's color will be changed into white. Otherwise, black. After this process, the image turns out to be dichromatic image with only black and white color. Although the process may lose some information, it is an indispensable step of some processes.




The curve on the **Binary** dialog shows different gray levels of the image. The line in the dialog indicates the threshold value. To change the value, you can drag it to change the threshold, or change the number in the Edit control on the dialog.

Click the **Best Fit** button to apply the auto threshold process to the image. The process offers an automatic threshold, click **Best Fit** can help make the image **Binary**.

### 5.5 Emboss... Shift+E

Operating method: Choose **Process>Emboss** menu





**Emboss** is a kind of artistic process, which makes the image look like an embossed image. The **Preview** button in **Emboss** dialog allows previewing the image before creating it. The process supplies 3 kinds of convolutions including **Gradient**, **Different**, and **Prewitt**. There are 8 directions for each convolution method. Users can get different effects with different convolution methods or directions.

## 5.6 Pseudo Color...

Operating method: **Process>Pseudo Color** menu.





Image must be in **Gray** mode. **Pseudo Color** command colorizes the active monochromatic image so as to highlight certain features in a gray scale image. The command does not modify the pixels' values in image bitmap (it does not convert image to true color or palette.), but creates a special palette with which the monochromatic image is displayed. It simply associates a **Pseudo Color** palette with the image that interprets the gray-level values in the image as color.

To select the start and end colors of the range, please click the buttons at each end of the color strip separately to bring up the color dialogues. The set result displays at the middle gradient bar.

## 5.7 Surface Plot...

Operating method: Choose **Progress>Surface Plot...** menu.




This command tool reflects the image intensity values intuitively as a 3-D, while in 3-D coordinate, X axes represents length, Y axes represents width, Z axes represents brightness. The rendering effect can be achieved in the preview window, and image position can be drawn in the window by adjusting the mouse.

Item	Description
3-D Preview Window Direction	<ul style="list-style-type: none"> <li>◆ Setting the edit box on the left side can adjust the 3-D surface plot position in the preview window, the default value is -0.5</li> <li>◆ Setting the edit box on the right side can adjust the relative. It's default value is 1</li> </ul>
Image Background Color	◆ Adjust the preview window <b>Image Background Color</b> by activating the <b>Color</b> dialog box system
Capture	◆ <b>Capture</b> the active image in the viewpoint window as a new image
Color Table	◆ Choose an appropriate color to display the different gray values in this table. Click the left color button to set the start color value, click the right color to set the end color value, the intermediate color bar displays the changes between the start point color and the end point( the relative information can be referred as color)
Reset	◆ Set the <b>Position Solid</b> 's two edit controls to their default values

## 5.8 Line Profile...

Operating method: Choose **Process>Line Profile** menu.




This tool illustrates how pixels along a selected line are distributed by graphing the number of pixels at each color intensity level. In a Line Profile, the X-axis represents the spatial scale, and the Y-axis represents the intensity values which range from 0 to 255.

Item	Description
Background	◆ Invoke the windows <b>Color</b> dialog to set the <b>Background Color</b> of the profile window.
Title	◆ Use this command to set a <b>Title</b> on the <b>Line Profile</b> image

Capture	◆ Capture the image in the <b>Line Profile</b> window as a new untitled image.
Copy	◆ Copy the <b>Line Profile</b> window's content onto the clipboard.
Save as	◆ Save the <b>Line Profile</b> image in bmp format.

## 5.9 Diffuse•••Shift+D

Operating method: Choose **Process>Diffuse** menu.




**Diffuse** is a kind of artistic process. It can diffuse the image. Users can adjust the parameter in the dialog to control the degree of the diffusion. Value range: odd number between 1~30.

Item	Description
Size	◆ Illustrate the diffuse degree, default value:1,Range:1~29
Preview	◆ Check it to display the real-time effect when drag the slider bar

## 5.10 Granulate•••Shift+G

Operating method: Choose **Process>Granulate** menu



**Granulate** is a process that can make the image blur. User can change the degree of image blurring, the value range:1~20.

Item	Description
Size	◆ Illustrate the <b>Granulate</b> size, default value:1,Range:1~20
Preview	◆ Check it to display the real-time effect when drag the slider bar

## 6 Copyright & Its Term for USage

### 6.1 Copyright Statement

**TOUPTEK PHOTONICS** is **TouView** software copyright owner, developer and owner, its software CD and accompanying documents are protected by international copyright laws. Our company will pursue its legal responsibility for illegal use and behavior in violation of copyright law. **TOUPTEK PHOTONICS** allow you backup software to protect your investment from loss, but only as a file for backup purposes, cannot be used for other commercial uses and purposes. You can not in any form of my company's software to re-distribution, modification, distribution, translate, decompile, tampering, crack and attempt to discover the source code of the software. **ACCU-View** cannot be provided for the third parties any writing or making and using of very similar software. Sequences license number is the only proof of your identification on the CD sleeve.

### 6.2 Terms for usage

#### **TOUPTEK SOFTWARE LICENSE TERMS**

These license terms are an agreement between **ToupTek** Corporation (or based on where you live, one of its affiliates) and you. Please read them. They apply to the software named above, which includes the media on which you received it, if any. The terms also apply to any **ACCU-View**

- updates,
- supplements
- Internet-based services, and
- support services

for this software, unless other terms accompany those items. If so, those terms apply.

**BY USING THE SOFTWARE, YOU ACCEPT THESE TERMS. IF YOU DO NOT ACCEPT THEM, DO NOT USE THE SOFTWARE.**

If you comply with these license terms, you have the rights below.

1. **INSTALLATION AND USE RIGHTS.**

- **Installation and Use.** You may install and use any number of copies of the software on your devices.
- **Included Microsoft Programs.** The software contains other Microsoft programs. The license terms with those programs apply to your use of them.

2. **Scope of License.** The software is licensed, not sold. This agreement only gives you some rights to use the software. ToupTek reserves all other rights. Unless applicable law gives you more rights despite this limitation, you may use the software only as expressly permitted in this agreement. In doing so, you must comply with any technical limitations in the software that only allow you to use it in certain ways. You may not

- disclose the results of any benchmark tests of the software to any third party without ToupTek's prior written approval;
  - work around any technical limitations in the software;
  - reverse engineer, decompile or disassemble the software, except and only to the extent that applicable law expressly permits, despite this limitation;
  - make more copies of the software than specified in this agreement or allowed by applicable law, despite this limitation;
  - publish the software for others to copy;
  - rent, lease or lend the software; or
  - use the software for commercial software hosting services.
3. **BACKUP COPY.** You may make unlimited number of backup copies of the software. You may use it only to reinstall the software.
4. **DOCUMENTATION.** Any person that has valid access to your computer or internal network may copy and use the documentation for your internal, reference purposes.
5. **TRANSFER TO A THIRD PARTY.** The first user of the software may transfer it and this agreement directly to a third party. Before the transfer, that party must agree that this agreement applies to the transfer and use of the software.
6. **USE WITH A THIRD PARTY CAMERA.** The use of this software with a third party camera is strictly prohibited. Any problems caused by use of this software with a third party camera are not supported by ToupTek. ToupTek reserves the right to investigate the behavior of using the software with third party cameras.
7. **EXPORT RESTRICTIONS.** The software is subject to P.R. China export laws and regulations. You must comply with all domestic and international export laws and regulations that apply to the software. These laws include restrictions on destinations, end users and end use.
8. **SUPPORT SERVICES.** Because this software is as is, we may not provide support services for it.
9. **ENTIRE AGREEMENT.** This agreement, and the terms for supplements, updates, Internet-based services and support services that you use, are the entire agreement for the software and support services.
10. **APPLICABLE LAW.**
- P.R. China. If you acquired the software in the United States, Washington state law governs the interpretation of this agreement and applies to claims for breach of it, regardless of conflict of laws principles. The laws of the state where you live govern all other claims, including claims under state consumer protection laws, unfair competition laws, and in tort.
  - Outside the P.R. China. If you acquired the software in any other country, the laws of that country apply.
11. **LEGAL EFFECT.** This agreement describes certain legal rights. You may have other rights under the laws of your country. You may also have rights with respect to the party from whom you acquired the software. This agreement

does not change your rights under the laws of your country if the laws of your country do not permit it to do so.

12. **DISCLAIMER OF WARRANTY.** The software is licensed as-is. You bear the risk of using it. ToupTEK gives no express warranties, guarantees or conditions. You may have additional consumer rights under your local laws which this agreement cannot change. To the extent permitted under your local laws, touptek excludes the implied warranties of merchantability, fitness for a particular purpose and non-infringement.
13. **LIMITATION ON AND EXCLUSION OF REMEDIES AND DAMAGES.** You can recover from touptek and its suppliers only direct damages up to U.S. \$5.00. You cannot recover any other damages, including consequential, lost profits, special, indirect or incidental damages.

This limitation applies to

- anything related to the software, services, content (including code) on third party Internet sites or third party programs;
- An claims for breach of contract, breach of warranty, guarantee or condition, strict liability, negligence, or other tort to the extent permitted by applicable law.

It also applies even if ToupTek knew or should have known about the possibility of the damages. The above limitation or exclusion may not apply to you because your country may not allow the exclusion or limitation of incidental, consequential or other damages.