



Bulb Replacement Guide

Some of our microscope models are designed with specific stands, while others can be configured on any number of our stands.

Be sure to select the bulb for your specific microscope or stand.

If you have any questions, contact your authorized ACCU-SCOPE dealer, or call us at 631-864-1000.

Model Series	Bulb Catalog #	Bulb Specifications
3000	#3256	6 volt 20 watt halogen
3000-LED	#00-3258-LED	3 watt LED
3002	#3256	6 volt 20 watt halogen
	#02-3258	LED light assembly
3002-CB	#3256	6 volt 20 watt halogen
	#02-3258	LED light assembly
3004	#3256	6 volt 20 watt halogen
3012	#3253	12 volt 20 watt halogen
	#12-3258-LED	3 watt LED
3015	#3256	6 volt 20 watt halogen
3025	#3254	Osram HBO 100W/2 mercury
	#3256-30W	6 volt 30 watt halogen
	#25-3252-G4	24 volt 100 watt halogen for dual view models
	#25-3252-G6	24 volt 100 watt halogen for Item #25-3232
3030	#3256-35	6 volt 20 watt halogen
3031	#3254	Osram HBO 100W/2 mercury
	#3256-35	6 volt 20 watt halogen
	#3256-30W	6 volt 30 watt halogen
3032	#3254	Osram HBO 100W/2 mercury
	#3256-30W	6 volt 30 watt halogen
3035	#3256-35	6 volt 20 watt halogen
3050	#3367	12 volt 10 watt fuse-style
3055	#3363	12 volt 10 watt halogen
3061 (top)	#3370	6 volt 15 watt dichroic halogen
3061 (top S/N 99 only)	#3369	6 volt 10 watt halogen
3061 (bottom)	#3368	5 watt fluorescent
3065 (top)	#3370	6 volt 15 watt dichroic halogen
3065 (bottom)	#3368	5 watt fluorescent
3069 (top)	#3370	6 volt 15 watt dichroic halogen
3069 (pre-2000 only)	#3369	6 volt 10 watt halogen
3069 (bottom)	#3368	5 watt fluorescent
3075	#3373	8 watt fluorescent bulb for Item #3360
	#3378	EKE 21 volt 150 watt, MR-16 lamp for Item #16105, 16106, or 16108 illuminators
3079	#3368	5 watt fluorescent
3079 (S/N 95-99 only)	#3276	20 watt oblong clear round
3079 (S/N 00 only)	#3260-25	25 watt screw
3080	#3257	7 watt fluorescent
3088F	#3368	5 watt fluorescent
	#88-3258-LED	LED



Bulb Replacement Guide

Model Series	Bulb Catalog #	Bulb Specifications
3359 Nicholas Style Illuminator	#3372	12 volt 15 watt dichroic halogen
3374 Fiber Optic Illuminator	#3378	20 volt 150 watt DDL lamp with reflector
3360 Fluorescent Ring Illuminator.....	#3373	8 watt fluorescent ring
3079 (S/N 95-99 only)	#3276	20 watt oblong clear round
3079 (S/N 00 only).....	#3260-25	25 watt screw
3085 (S/N 95-99 only)	#3276	20 watt oblong clear round
3085 (post 2000 only).....	#3257-85F	7 watt fluorescent