



**3050 Stereo Microscope Series**

<b>Catalog No.</b>	<b>Description</b>
<b>3050-LED</b>	<b>Stereo Microscope with 1x and 2x objectives</b> WF10x/20mm eyepieces Frosted glass and black & white stage plates, eyeguards & dust cover Rechargeable LED reflected and transmitted illumination
<b>3051-LED</b>	<b>Stereo Microscope with 2x and 4x objectives</b> WF10x/20mm eyepieces Frosted glass and black & white stage plates, eyeguards & dust cover Rechargeable LED reflected and transmitted illumination
<b>3052-LED</b>	<b>Stereo Microscope with 1x and 3x objectives</b> WF10x/20mm eyepieces Frosted glass and black & white stage plates, eyeguards & dust cover Rechargeable LED reflected and transmitted illumination
<b>3053-LED</b>	<b>Stereo Microscope with 1x and 4x objectives</b> WF10x/20mm eyepieces Frosted glass and black & white stage plates, eyeguards & dust cover Rechargeable LED reflected and transmitted illumination

<b>Eyepieces</b>	
3306-50	WF10x stereo eyepiece, 20mm field of view, 30.5mm O.D. (each)
3307-50	WF15x stereo eyepiece, 13mm field of view, 30.5mm O.D. (each)
3308-50	WF20x stereo eyepiece, 10mm field of view, 30.5mm O.D. (each)
3319	Stereo eyeguards (pair)

<b>Miscellaneous Accessories</b>	
1111	Microscope carry case, hard-sided with fabric interior, accessory pockets, and velcro straps Interior Dimensions: 10.5" L x 7.5"W x 13.5"H
3301-S	Dust cover
3346-50	Darkfield condenser
3347	Wire gem clamp with mounting post
3350-50	Frosted glass stage plate for 3050, 95mm diameter
3351-50	Black and white stage plate for 3050, 95mm diameter
3367	12 volt 10 watt fuse style bulb, 5/pk