

3075 ZOOM STEREO MICROSCOPE SERIES

QUICK REFERENCE GUIDE

The ACCU-SCOPE 3075 Zoom Stereo Microscope Series provides sharp, high resolution, three-dimensional images in a powerful professional microscope designed to withstand rigorous industrial use. The large magnification range and long working distance simplifies observation making it ideal for industrial inspection applications.

Viewing Heads



45° Inclined
Binocular Head
#75-3330



45° Inclined
Trinocular Head
#75-3332

Viewing Heads

Item#	Description
75-3330	Inclined 45° Binocular
75-3332	Inclined 45° Trinocular

Eyepieces



WF10x/22mm FOV
#75-3311 & #75-3311R



WF15x/16mm FOV
#75-3312



WF20x/12mm FOV
#75-3313

Eyepieces

Item#	Description
75-3311	10x/22mm Eyepiece
75-3311-R	10x/22mm Eyepiece (With 10/100mm reticle)
75-3312	15x/16mm Eyepiece
75-3313	20x/12mm Eyepiece

Objectives - Auxiliary



0.3x/WD 287mm
#75-3333



0.4x/WD 217mm
#75-3334



0.5x/WD 177mm
#75-3336



0.7x/WD 120mm
#75-3337



1.0x/WD 100mm
Glass Lens Shield
#75-3338



1.5x/WD 47mm
#75-3339



2.0x/WD 26mm
#75-3340

Auxiliary Objectives

Item#	Description
75-3333	0.3x
75-3334	0.4x
75-3336	0.5x
75-3337	0.7x
75-3338	1.0x Glass Lens Shield
75-3339	1.5x
75-3340	2.0x

Adapters



0.35x C-mount
#75-2010-35



0.50x C-mount
#75-2010-50



0.65x C-mount
#75-2010-65



Photo Tube - 23.2mm ID
#75-2010-P

Adapters

Item#	Description
75-2010-35	0.35x C-mount Adapter
75-2010-05	0.50x C-mount Adapter
75-2010-65	0.65x C-mount Adapter
75-2010-P	Photo Tube

V041718