

3078-HDR Zoom Stereo Microscope



3078-HDR Zoom Stereo Microscope

Catalog No.	Description
3078-HDR	<p>Binocular Zoom Stereo Microscope With Built-In HD Camera Binocular viewing head, inclined 60° WF10x/20mm focusable eyepieces with eyeguards 0.8x to 3.5x zoom range (8x -35x with 10x eyepieces) 5 mega pixel built-in CMOS camera 1080p HD output and AV output, 4 GB SD card for image capture USB 2.0 PC compatible (<i>image transfer only - no live preview on a computer</i>) Remote control, HDMI, USB, AV cables all included Dovetail focusing stand with coarse focus adjustment & stage plate Variable 5 watt LED incident and transmitted illuminators, 100v - 240v universal power supply</p>

Eyepieces	
75-3311	WF10x/22mm eyepieces with built in diopter adjustment & eyeguard
75-3311-R	WF10x/22 focusable eyepiece with 10mm/100 division reticle with cross-line
75-3312	WF15x/16mm eyepiece with built in diopter adjustment & eyeguard
75-3313	WF20x/12mm eyepiece with built-in diopter adjustment & eyeguard
75-3319	Stereo eyeguards (pair)