



3080 Microscope Series

Catalog No.	Description
3080	Monocular Microscope with Disc Diaphragm Monocular viewing head, inclined 45°, rotatable 360° WF10x/18mm eyepiece with pointer Reversed quadruple nosepiece DIN 4x, 10x and 40xR objectives Flat top stage with stage clips, disc diaphragm Separate coarse and fine focusing controls 7 watt fluorescent illuminator and cord hanger

Viewing Heads, Eyepieces and Reticles

3150-80	Monocular viewing head, inclined 45°
3105-80	WF10x/18mm field of view eyepiece with pointer
3106	WF15x/15mm field of view eyepiece
3107	WF20x/12mm field of view eyepiece

Objectives

3195-80	4x achromat objective, N.A. 0.10, W.D. 35.8mm
3196-80	10x achromat objective, N.A. 0.25, W.D. 6.47mm
3197-00	20x achromat objective, N.A. 0.40, W.D. 1.0mm
3198-80	40xR achromat objective, N.A. 0.65, W.D. 0.45mm
02-3199	60xR DIN achromat objective, N.A. 0.85, W.D. 0.19mm

Optional Accessories

3214-80	Stage clips, locked-on w/spring, M3 x 14, 55mm blade length
3257	7 watt fluorescent bulb, 3/pk
3301-S	Dust cover - 19.5" H x 15.5" W
3277-80	0.8 amp fuse, 5/pk