

Catalog No.	Description
Compound Eyepieces & Roll Down Eyeguards (23.2mm I.D.)	
3104	WF10x eyepiece, 18mm field of view
3105	WF10x eyepiece with pointer, 18mm field of view
02-3104	WF10x/18.5mm field of view eyepiece with reticle holder (3002 Series)
02-3105	WF10x/18.5mm field of view eyepiece with pointer & reticle holder (3002 Series)
3110-00	HWF10x eyepiece, 18mm field of view (3000 Series S/N 02 - 11xxxx)
3111-00	HWF10x eyepiece with pointer, 18mm field of view (3000 Series S/N 01-11xxxx)
3110-15	HWF10x eyepiece, 20mm field of view (3015 Series)
3106	WF15x eyepiece, 15mm field of view
3107	WF20x eyepiece, 11mm field of view
3131	Eyeguards, roll down (pair) for 23.2mm I.D. eyepiece

Compound Eyepieces & Roll Down Eyeguards (30mm I.D.)	
12-3110	HWF10x eyepiece, 22mm field of view with built in diopter adjustment. Eyepiece accepts a 27mm reticle and has a 30mm O.D. (3012 Series)
25-3113	EWf 10x/22mm focusable eyepiece with diopter adjustment (3025 Series)
25-3113-CL	EWf 10x/22mm focusable eyepiece with 10mm/100 division crossline reticle
25-3114	EWf 10x/20mm focusable eyepiece with diopter adjustment
12-3106	WF15x eyepiece with built-in diopter adjustment, 16mm field of view. Eyepiece accepts a 20mm reticle and has a 30mm O.D.
12-3107	WF20x eyepiece with built-in diopter adjustment, 12.5mm field of view. Eyepiece does not accept a reticle and has a 30mm O.D.
25-3131	Roll down eyeguards (pair)

Pointers, Reticles & Stage Micrometers	
3130	Pointers, spring-wire, 50/pk. - for 23.2mm I.D. eyepieces
12-3130	27mm Eyepiece reticle pointer (each), for 30mm I.D. eyetubes
3137-19	19mm eyepiece reticle w/mounting ring, 10mm/100 divisions with crossline
3137-19-PS	19mm eyepiece reticle with pointer and scale
3137-21	21mm eyepiece reticle w/mounting ring, 10mm/100 divisions with crossline
25-3137-26	26mm eyepiece reticle w/mounting ring, 10mm/100 divisions with crossline
3145	Stage micrometer, 1mm/100 divisions (each division .010mm)

DIN Plan Achromat Objectives	
3193SP-00	2x DIN semi-plan achromat objective, N.A. 0.05, W.D. 17.0mm
3185PL-00	4x DIN plan achromat objective, N.A. 0.10, W.D. 17.9mm
3186PL-00	10x DIN plan achromat objective, N.A. 0.25, W.D. 2.04mm
3187PL-00	20xR DIN plan achromat objective, N.A. 0.40, W.D. 1.2mm
3188PL-00	40xR DIN plan achromat objective, N.A. 0.65, W.D. 0.65mm
3190PL-00	50xR DIN oil plan achromat objective, N.A. 0.95, W.D. 0.18mm
3191PL-00	60xR DIN plan achromat objective, N.A. 0.85, W.D. 0.19mm
3189PL-00	100xR DIN plan achromat objective, N.A. 1.25, W.D. 0.15mm

Catalog No.	Description
Non-DIN Objectives for Bausch & Lomb Galen II® & Bristoline®	
3165	4x Non-DIN achromat objective, N.A. 0.10, W.D. 34.70mm
3166	10x Non-DIN achromat objective, N.A. 0.25, W.D. 7.63mm
3168	40xR Non-DIN achromat objective, N.A. 0.65, W.D. 0.53mm
3170	100xR Non-DIN achromat objective, N.A. 1.25, W.D. 0.198mm
3166A	A10x Non-DIN semi-plan objective, N.A. 0.10, W.D. 6.6mm (Japanese)

Catalog No.	Description
Objectives for the 3017 Series	
3195SP-17	4x DIN semi-plan objective for 3017, N.A. 0.10, W.D. 35.8mm
3196SP-17	10x DIN semi-plan objective for 3017, N.A. 0.25, W.D. 2.49mm
3198SP-17	40xR DIN semi-plan objective for 3017, N.A. 0.65, W.D. 0.65
3200SP-17	100xR DIN semi-plan oil objective for 3017, N.A. 1.25, W.D. 0.52mm

Catalog No.	Description
Objectives for the 3085 Series	
3196-85	10x DIN achromat objective for 3085, N.A. 0.25, W.D. 7.14mm
<i>All other objectives are discontinued. Use objectives from the 3088F Series as replacements</i>	

Catalog No.	Description
Parfocal Rings / Shims	
3206-1	.001" Shims - stainless steel, 10/pk.
3206-2	.002" Shims - stainless steel, 10/pk.
3206-3	.003" Shims - stainless steel, 10/pk.
3206-4	.004" Shims - stainless steel, 10/pk.
3206-5	.005" Shims - stainless steel, 10/pk.
3206-7	.007" Shims - stainless steel, 10/pk.
3206-10	.010" Shims - stainless steel, 10/pk.

Catalog No.	Description
Attachable Mechanical Stages	
3210	Attachable mechanical stage - side controls, straight mount (Chinese)
3210-J	Attachable mechanical stage - side controls, off-center mount (Japanese style)
3211	Attachable mechanical stage - top controls - straight mount (Chinese)
3211-J	Attachable mechanical stage - top controls, off center mount (Japanese style)

Catalog No.	Description
Stage Clips and Finger Springs	
3214	Stage clips, push-in, M3 x 5, 40mm blade length (pair)
3214-80	Stage clips, locked-on w/spring, M3 x 14, 55mm blade length (pair)
3214-88	Stage clip, locked-on w/spring, M3 x 8, blade length 55mm (pair)
3215-00	Finger spring for 3000 Microscope Series
3215-04	Finger spring for 3004 Microscope Series
3215-15	Finger spring for 3015 Microscope Series
02-3215	Finger spring for 3002 Microscope Series
12-3215	Finger spring for 3012 Microscope Series
3216-J	Finger spring, left for Bristol®/Japanese
3217-J	Finger spring, right for Bristol®/Japanese
3218	Finger spring/attachable mechanical stage

Catalog No.	Description
Double Plate Mechanical Stages and Slide Holders	
02-3213	Double plate mechanical stage for 3002 Series
00-3213	Double plate mechanical stage for 3000-LED series
12-3213	Double plate, rackless mechanical stage for the 3012 Series
3213-00	Double plate mechanical stage for 3000 Series
3213-00-200	Double plate mechanical stage for 3000 Series (S/N 94 thru 01)
3213-04	Double plate mechanical stage for 3004 Series
3213-15	Double plate mechanical stage for 3015 Series
3213-88	Double plate mechanical stage for 3088F Series
02-SH	Slide holder for 3002 Series mechanical stage
00-SH	Slide holder for 3000-LED Series mechanical stage
12-SH	Slide holder for 3012 Series mechanical stage
SH-00	Slide holder for 3000 Series mechanical stage
SH-00-200	Slide holder for 3000 Series mechanical stage (S/N 94 thru 01)
SH-04	Slide holder for 3004 Series mechanical stage
SH-15	Slide holder for 3015 Series mechanical stage
SH-88	Slide holder for 3088F Series mechanical stage

Condensers	
3221-J	N.A. 1.25 Abbe condenser w/iris diaphragm & filter holder for Bristoline® / Japanese microscope, Ø 36mm

Filters	
25-3240	45mm Clear daylight blue filter
32-3240-G	45mm Green filter
3235	32mm Clear daylight blue filter
3235-G	32mm Green filter
3235-W	32mm White frosted filter
3236	32mm Frosted daylight blue filter
3237	37mm Clear daylight blue filter
3238	37mm Frosted daylight blue filter
3238-15	40mm Clear daylight blue filter
12-3222	Filter holder with blue filter for 3012 Series <i>also works on Leica® CME</i>

Attachable Illuminators, Converters & Cord Hangers	
3246	20 watt attachable illuminator with blue filter, front or bottom attachment, 5mm stem
3249	Mirror & fork (mirror dia. 50mm, post dia. - 5mm, post length - 12mm)
3250	120/220v Converter
3251	Cord hanger with screws (pair)

Catalog No.	Description
Bulbs	
3253	Osram® #64425, 12 volt 20 watt halogen - 5/pk.
3254-U	Ushio® HBO 100 watt mercury bulb
3255	6 volt 3 amp bayonet bulb, 5/pk.
3256	Osram® #64250, 6 volt 20 watt halogen bulb - 5/pk.
3256-30W	Philips® #5761, 6 volt 30 watt halogen bulb - 5/pk.
3257	7 watt fluorescent bulb, 3/pk., (3080, 3085F post 2000)
3368	5 watt fluorescent bulb (2 pin style), 3/pk. (3088F, pre-2000 3085F)
3260-20	20 watt screw bulb, 5/pk.
3260-25	25 watt screw bulb, 5/pk. (3079 S/N 00)
3260-30	30 watt screw bulb, 5/pk.
3262	30 watt bayonet mirror bulb, 5/pk. (Galen II® and Bristol®)
3268	12 volt 2 amp bayonet bulb, 5/pk.
3275	20 watt bayonet bulb clear round, 5/pk.
3276	20 watt oblong bulb clear round, 5/pk. (3085, 3079 S/N 95 - 99)

Circuit Boards and Power Supplies	
3280	4-wire capacitor for Bristol® 3000
3281	5-wire capacitor for Bristol® 3000
3282	American Optical® Series 10 lamp socket (use w/A.O. 6v3a bulb# 1460)
3283	American Optical® Series 10 bulb installer for A.O. bulb #1460 (6v3a)
3284	Circuit board direct wire for 3000
3293-04	Replacement power supply with rheostat for 3004 Series
3291	Circuit board w/male adapter for Bristol® 3002
3292	Circuit board direct wire for Bristol® 3002
3293	Black box - complete replacement power supply w/ wiring diagram. Includes transformer, circuit board, socket, rheostat & 6v20w halogen bulb
3293-30	30 watt external power supply with variable intensity control (use for replacements on older microscopes)

Sockets and Rheostats	
3290-SQ	Lamp socket - halogen, square (for 6 volt 20 watt halogen bulb)
3290-R	Lamp socket - halogen, round (for 6 volt 20 watt halogen bulb)
3290-15	Lamp socket - halogen, square for 3015 series
02-3290	Lamp socket - square for 3002 Series
25-3290-30W	Lamp socket w/illuminator housing for 3025 series - 6v30w version
25-3290-50W	Lamp socket w/illuminator housing for 3025 series - 12v50w version
32-3290	Lamp socket w/illuminator housing for 3032 series - 6v30w version
G2-R	Rheostat for Bausch & Lomb Galen II® & Bristoline® microscopes
00-3287	Rheostat for 3000-LED Series
02-3287	Rheostat for 3002 Series
12-3287	Rheostat for 3012 Series
15-3287	Rheostat for 3015 Series
25-3287-UP	Rheostat for 3025/3032 Series - Universal Power Supply
25-3287	Rheostat for 3025 / 3032 Series - Switchable Power Supply

Catalog No.	Description
Microscope Carry Cases & Dust Covers	
1111	Microscope carry case, hard-sided with fabric interior, accessory pockets, and velcro straps Interior Dimensions: 10.5" L x 7.5"W x 13.5"H
3294	Microscope carry case, hard sided with fabric lining, accessory pockets and velcro straps. Interior dimen. 11"L x 7.75" W x 14.25" H
3295	Microscope carry case, hard side with fabric lining, accessory pockets and velcro straps. Interior dimen. 11.5" L x 9.25" W x 17.5" H
02-3295-A	Aluminium hard-sided carry case for 3002 Series (form fit)
00-3295-A	Aluminium hard sided carry case, form fitted for 3000/3001-LED
3296-15	Aluminum hard sided carry case, form fitted for the 3015 microscope
3301-S	Dust cover, black nylon - 15.5" W x 19.5" H
3301-L	Dust cover, black nylon - 18" W x 22" H
3301-XL	Dust cover, black nylon - extra large size
25-3301	Dust cover for 3025 Microscope
25-3301-2	Dust cover for 3025-2 Dual Viewing Microscope
25-3301-2F	Dust cover for 3025-2F Front/Back Dual Viewing Microscope
25-3301-5	Dust cover for 3025-5 Multi-Viewing Microscope

Prepared Slide Sets - Biology	
3400	15 piece assorted prepared slide set in a box
3405	25 piece assorted prepared slide set in a box
3410	50 piece assorted prepared slide set in a box
3415	89 piece assorted prepared slide set in a box