

Catalog No.	Description
Bausch & Lomb® and Leica® Parts, Stands and Focus Mounts	
3115	WF10x eyepiece similar to B&L® # 31-15-81 for Stereozoom 4 and 6
3131	Eyeguards, roll-down (pair)
3338-BL	Glass lens shield for B&L® Stereozoom, similar to # 31-26-21
15350	Swing-arm stand similar to Leica® # 31-26-10, lockable 14.5" vertical motion range, 3 foot scanable range, adjustable counterbalanced support use with CAT# 3506-BL
15351	Horizontal or vertical clamp mount similar to Leica # 31-26-14, use w/Cat # 15350
15350-TM	Table mount for flexible arm stand (CAT# 15350); attaches w/4 screws (not included)
15350-WM	Wall mount for flexible arm stand (CAT# 15350); attaches w/4 screws not included
3354-BL	E-Arm (Bonder arm) similar to Bausch & Lomb® # 31-26-59
3354-L	E-Arm (Bonder arm) with 76mm yoke similar to Leica® # 10446345
3357	Universal (Boom) stand similar to Leica® S-Stand # 31-27-14
3357-BB	Ball bearing boom stand, horizontal arm with bonder port, base dimensions 10" x 10", Vertical post 15.5", Horizontal arm 23"

Stereo Stage Plates	
3350-50	Frosted glass stage plate for 3050, 95mm diameter
3350-55	Frosted glass stage plate, 59mm diameter
3350-61	Frosted glass stage plate for 3059/3061/3069 - 74mm diameter
3350-65	Frosted glass stage plate for 3065, 95mm diameter
3351-50	Black and white stage plate for 3050, 95mm diameter
3351-55	Black and white stage plate, 59mm diameter
3351-61	Black and white stage plate for 3059/3061/3069- 74mm diameter
3351-62	Black and white stage plate for 3062/3060, 89mm diameter
3351-65	Black and white stage plate for 3065, 95mm diameter
3353	Stage clip, spring loaded locked on, M4 x 6, 55mm blade length

Stereo Illuminator Bulbs	
3364	115 volt 15 watt screw stereo bulb, 5/pk.
3365	6 volt 1.2 amp top stereo bulb, 5/pk.
3366	6 volt 1.2 amp bottom stereo bulb, 5/pk.
3367	12 volt 10 watt fuse style bulb, 5/pk. (3050 Microscope Series)
3368-61	5 watt fluorescent bulb, 4 pin style (all models S/N 07+, bottom illum.)
3368	5 watt fluorescent bulb (2 pin style), 3/pk. (3061, 3069, 3065 Series), 3/pk.
3369	Osram® # 64225 - 6 volt 10 watt halogen bulb, 3/pk. (pre-2000 3061 & 3069 Series)
3370	6 volt 15 watt dichroic halogen, MR11 (3061, 3069, 3065 Series') - 3/pk.
3372	12 volt 15 watt dichroic halogen bulb, MR11 for Cat# 3359 - 3/pk.
3373	8 watt fluorescent ring illuminator bulb, 2/pk.