

3072 STEREO MICROSCOPE SERIES

QUICK REFERENCE GUIDE

Advanced, modular design incorporating many professional features the 3072 Stereo Microscope Series is ideal for classrooms, light industrial applications and hobbyists. A large working distance, enhanced illumination and eyepieces providing a large field of view make these microscopes a great choice for classroom dissections and industrial inspections.

Viewing Heads



1x/3x—45° Inclined
Binocular Head
#72-3330-13



1x/2x—45° Inclined
Binocular Head
#72-3330-12



2x/4x—45° Inclined
Binocular Head
#72-3330-24

Viewing Heads

Item#	Description
72-3330-13	Inclined 45° Binocular 1x/3x fixed objectives
72-3330-12	Inclined 45° Binocular 1x/2x fixed objectives
72-3330-24	Inclined 45° Binocular 2x/4x fixed objectives

Eyepieces



WF10x/20mm FOV
#132-10-11



WF15x/16mm FOV
#131-10-15



WF20x/12.5mm FOV
#131-10-20

Eyepieces

Item#	Description
132-10-11	10x/20mm Eyepiece
131-10-15	15x/16mm Eyepiece
131-10-20	20x/12.5mm Eyepiece

Objectives - Auxiliary



0.5x/Wd 180mm
Auxiliary
#132-15-05



0.63x/Wd 140mm
Auxiliary
#132-15-06



0.75x/Wd 120mm
Auxiliary
#132-15-75



2.0x/Wd 30mm
Auxiliary
#132-15-02

Auxiliary Objectives

Item#	Description
132-15-05	0.5x (Use with #131-20-40)
132-15-06	0.63x
132-15-75	0.75x
132-15-02	2.0x

Adapters



0.5x Auxiliary Lens
Riser Adapter
#131-20-40

Adapters

Item#	Description
131-20-40	Riser Adapter (Use with #132-15-05)