



3072 Stereo Microscope Series

Catalog No.	Description
3072-13	Binocular Stereo Microscope 1x/3x on E-LED Stand Binocular viewing head, inclined 45°, rotatable 360° WF10x/20mm eyepieces w/built in diopter adjustments & eyeguards 1x and 3x objective magnifications 1.0x glass lens shield Dovetail focusing stand with coarse focus adjustment Frosted glass & black & white stage plates Variable LED incident and transmitted illuminators
3072-12	Binocular Stereo Microscope 1x/2x on E-LED Stand Same as above except with 1x and 2x objective magnifications
3072-24	Binocular Stereo Microscope 2x/4x on E-LED Stand Same as above except with 2x and 4x objective magnifications
3072-13-PFS	Binocular Stereo Microscope 1x/3x on Plain Focusing Stand Binocular viewing head, inclined 45°, rotatable 360° WF10x/20mm eyepieces w/built in diopter adjustments & eyeguards 1x and 3x objective magnifications 1.0x glass lens shield Dovetail focusing stand with coarse focus adjustment, black & white stage plate and dust cover
3072-12-PFS	Binocular Stereo Microscope 1x/2x on Plain Focusing Stand Same as above except with 1x and 2x objective magnifications
3072-24-PFS	Binocular Stereo Microscope 2x/4x on Plain Focusing Stand Same as above except with 2x and 4x objective magnifications
3072-13-PS	Binocular Stereo Microscope 1x/3x on Pole Stand Binocular viewing head, inclined 45°, rotatable 360° WF10x/20mm eyepieces w/built in diopter adjustments & eyeguards 1x and 3x objective magnifications 1.0x glass lens shield Pole stand with coarse focus adjustment; focus range 50mm Black and white stage plate
3072-12-PS	Binocular Stereo Microscope 1x/2x on Pole Stand Same as above except with 1x and 2x objective magnifications
3072-24-PS	Binocular Stereo Microscope 2x/4x on Pole Stand Same as above except with 2x and 4x objective magnifications
3072-13-BS	Binocular Zoom Stereo Microscope 1x/3x on Boom Stand Binocular viewing head, inclined 45°, rotatable 360° WF10x/20mm eyepieces w/built in diopter adjustments & eyeguards 1x and 3x objective magnifications 1.0x glass lens shield Boom stand, bonder arm and dust cover
3072-12-BS	Binocular Zoom Stereo Microscope 1x/2x on Boom Stand Same as above except with 1x and 2x objective magnifications
3072-24-BS	Binocular Zoom Stereo Microscope 2x/4x on Boom Stand Same as above except with 2x and 4x objective magnifications



3072 Stereo Microscope Series

Catalog No.	Description
Viewing Heads	
72-3330-13	Binocular viewing head, 1x and 3x objectives, 100mm working distance
72-3330-12	Binocular viewing head, 1x and 2x objectives, 100mm working distance
72-3330-24	Binocular viewing head, 2x and 4x objectives, 100mm working distance
Eyepieces	
132-10-11	WF10x focusing eyepiece, 20mm field of view. Eyepiece accepts a 26.5mm reticle,
131-10-15	WF15x focusing eyepiece, 16mm field of view. Eyepiece accepts a 20mm reticle
131-10-20	WF20x focusing eyepiece, 12.5mm field of view. Eyepiece does not accept a reticle
131-10-01	Rolldown eyeguards (pair) for CAT# 132-10-11
Reticles	
132-12-01	26mm eyepiece reticle, 10mm/100 division with crossline
131-12-05	26mm eyepiece reticle, 10mm/100 divisions without crossline
131-12-10	20mm eyepiece reticle, 10mm/100 divisions with crossline
131-12-11	20mm eyepiece reticle, 10mm/100 divisions
131-12-20	Stage micrometer, 1mm/100 division (each division .010mm)
131-12-21	Calibration slide
Auxiliary Objectives	
132-15-05	0.5x auxiliary objective, 180mm working distance (riser adapter required on LED and Plain Focusing Stands)
131-20-40	Riser adapter for 0.5x auxiliary objective
132-15-06	0.63x auxiliary objective, 140mm working distance
132-15-75	0.75x auxiliary objective, 120mm working distance
131-15-00	1.0x auxiliary objective, 100mm working distance (<i>required for use of fiber optic ring light</i>)
132-15-02	2.0x auxiliary objective, 30mm working distance
Focus Mounts and Bonder Arms	
75-3354	range
131-14-18	Focus mount for pole stand with \varnothing 32mm opening for post
Stages & Stage Plates	
131-20-05	Black and white stage plate, 125mm for Plain Focusing & Pole Stand
131-20-07	Black and white stage plate, 95mm for E-LED
131-20-14	Clear glass stage plate, 95mm for E-LED stand
131-20-21	Stage clips (pair) for plain focusing stand (CAT# 131-13-05)
3346-65	Darkfield condenser with adjustable iris diaphragm, 95mm O.D.
131-20-35	Wire gem clamp with mounting post
Miscellaneous	
3301-L	Dust cover