

## **ScopeLED® Fluorescence & Brightfield Illuminators**

---

### **ScopeLED® F-Series Fluorescence Illuminators**

CATALOG #	DESCRIPTION
SLF110-LX	ScopeLED single wavelength LED excitation fluorescence light source, control pod and power source - See attached ScopeLED wavelength chart to select wavelengths
SLF120-PX	ScopeLED F120 package - choice of two fixed wavelengths, control pod and power source - See attached ScopeLED wavelength chart to select wavelengths
SLF140-GX	ScopeLED F140 package - four fixed wavelengths, control pod, power source - See attached ScopeLED wavelength chart to select wavelengths

### **ScopeLED® Wavelength Options**

390±5, 440±5, 470±10, 525, 540±5, 560P(±10), 595±5, 640±5, 660±5, 740±5

### **ScopeLED® F-Series Adapters**

SLAF-LC1	L1 F-Series Leica Fluor Illuminator
SLAF-LT3	L3 F-Series Leitz Fluor Illuminator Old Style with Locking Ring
SLAF-NK1	N1 F-Series Nikon Std Fluor Illuminator Bayonette Mount
SLAF-NK2	N2 F-Series Nikon Fluor Illuminator TS100 E200
SLAF-OL1	O1 F-Series EXC-350/500 & Olympus Std Fluor Illuminator
SLAF-OL2	O2 F-Series Olympus BH Fluor Illuminator
SLAF-ZS1	Z1 F-Series Adapter for Zeiss Std Dovetail
SLAF-W1	W1 F-Series Adapter for VWR Vista Vision
SLAF-2BNC	F-Series 4 pin to 2 BNC Interface Cable Single Channel TTL/Analog