



EXS-210 Stereo Microscope Series

CATALOG #	DESCRIPTION
EXS-210-12	<p>Stereo Microscope with 1x and 2x objectives Binocular viewing head, inclined 45°, rotatable 360° WF10x/20mm eyepieces with eyeguards 1x and 2x objective magnifications Frosted glass and black & white plates, dust cover Variable reflected and transmitted rechargeable LED illumination for both corded and cordless operation</p>
EXS-210-24	<p>Stereo Microscope with 2x and 4x objectives Binocular viewing head, inclined 45°, rotatable 360° WF10x/20mm eyepieces with eyeguards 2x and 4x objective magnifications Frosted glass and black & white plates, dust cover Variable reflected and transmitted rechargeable LED illumination for both corded and cordless operation</p>
EXS-210-13	<p>Stereo Microscope with 1x and 3x objectives Binocular viewing head, inclined 45°, rotatable 360° WF10x/20mm eyepieces with eyeguards 1x and 3x objective magnifications Frosted glass and black & white plates, dust cover Variable reflected and transmitted rechargeable LED illumination for both corded and cordless operation</p>
EXS-210-123	<p>Stereo Microscope with 1x, 2x and 3x objectives Binocular viewing head, inclined 45°, rotatable 360° WF10x/20mm eyepieces with eyeguards 1x, 2x and 3x objective magnifications Frosted glass and black & white plates, dust cover Variable reflected and transmitted rechargeable LED illumination for both corded and cordless operation</p>
EXS-210-124	<p>Stereo Microscope with 1x, 2x and 4x objectives Binocular viewing head, inclined 45°, rotatable 360° WF10x/20mm eyepieces with eyeguards 1x, 2x and 3x objective magnifications Frosted glass and black & white plates, dust cover Variable reflected and transmitted rechargeable LED illumination for both corded and cordless operation</p>



EXS-210 Stereo Microscope Series

CATALOG #	DESCRIPTION	LIST PRICE
-----------	-------------	------------

Eyepieces

3306-50	WF10x stereo eyepiece, 20mm field of view, 30.5mm O.D. (each)
3307-50	WF15x stereo eyepiece, 13mm field of view, 30.5mm O.D. (each)
3308-50	WF20x stereo eyepiece, 10mm field of view, 30.5mm O.D. (each)
210-3319	Stereo eyeguards (pair)

Miscellaneous Accessories

1111	Microscope carry case, hard-sided with fabric interior, accessory pockets, and velcro straps Interior Dimensions: 10.5" L x 7.5"W x 13.5"H
3301-S	Dust cover
3350-55	Frosted glass stage plate for 3055, 59mm diameter
3351-55	Black and white stage plate for 3055, 59mm diameter
210-3365-B	3.5V 20mA LED replacement board set
210-3365-T	3.5 volt 1 watt LED replacement lamp
210-3259	100V-240V 4.5V 1A AC/DC Adapter