



Bulb Replacement Guide

Some of our microscope models are designed with specific stands, while others can be configured on any number of our stands. Be sure to select the bulb for your specific microscope or stand.

If you have any questions, contact your authorized ACCU-SCOPE dealer, or call 631-864-1000.

Model Series	Bulb Catalog #	Bulb Specifications
131-13-28 (<i>Diascopic stand</i>).....	131-13-28-LED	Replacement LED lamphouse
3000	3256.....	6 volt 20 watt halogen
3000-LED	00-3365	3 watt led
3002	3256.....	6 volt 20 watt halogen
	02-3258-SLED.....	LED light assembly
	02-SLED-KIT	Complete base, converts from halogen to LED
3002-CB.....	3256.....	6 volt 20 watt halogen
	02-3258-SLED.....	LED light assembly
3004	3256.....	6 volt 20 watt halogen
3012	3256.....	6 volt 20 watt halogen
	12-3258-LED	3 watt led
	12-3259-LED	5 watt led
3015.....	3256.....	6 volt 20 watt halogen
3025.....	3254-U.....	Ushio HBO 100 W/2 mercury
	25-3365	3 watt led
	25-3253-50	12 volt 50 watt
3031	3254-U.....	Ushio HBO 100 W/2 mercury
3032	3254-U.....	Ushio HBO 100 W/2 mercury
3052 (<i>REFLECTED</i>)	50-3258-LED-R	LED lamp and board
3052 (<i>TRANSMITTED</i>)	50-3258-LED-T.....	LED lamp and board
3061 (<i>TOP - Pre-2000 only</i>).....	3369.....	6 volt 10 watt halogen
3069 (<i>TOP - Pre-2000 only</i>).....	3369.....	6 volt 10 watt halogen
3078-HD	78-3365	5 watt led
3088-LED	88-3258-LED	Replacement LED
EXC-100	100-3258	2 watt replacement LED
EXC-120	120-3258-SLED.....	1 watt LED replacement
EXC-350	350-3365	3 watt led
EXC-360	360-3365	3 watt LED
EXC-400	400-3365	3 watt led
EXC-500	500-3290	Replacement lamphouse with LED bulb
EXI-310.....	310-3365	5 watt LED
EXI-410.....	410-3365	3 watt SLED replacement with illuminator housing
EXI-600.....	600-3290	10 watt LED light module
EXM-50.....	50-3258	Replacement LED
EXM-150.....	150-3258-LED	Replacement lamphouse with LED bulb
EXS-210 (<i>TOP</i>).....	210-3365-T	1 watt LED replacement
EXS-210 (<i>BOTTOM</i>).....	210-3365-B.....	Replacement LED