

See the TRUE and ACCURATE Representation of Your Samples!

Microorganisms often have very little to no coloration, so they are difficult to see with traditional microscopy methods.

Darkfield works differently. Any object in the sample, including microorganisms, will scatter light. This light enters the microscope objective and is seen as a bright object on a dark background (darkfield). The background is dark because there is nothing present to scatter the light into the objective.

View Your Samples in a Different Light!

Darkfield is available as an easy-to-use swing-in/swing-out accessory that mounts to the bottom of the condenser of your microscope.

The ACCU-SCOPE darkfield annulus (CAT# 150-3220) is available on the EXC-120 and EXC-100 series microscopes and select models of the EXM-150 series.

